

| Main Differences | |
|---------------------------|--|
| <p>Static SPR</p> | <p>Purpose: Scouting or screening prior to Kinetic SPR, yes/no binding</p> <ul style="list-style-type: none"> • Like label-free ELISA: ligand-analyte interaction depends on diffusion, concentration, and binding strength is measured as Apparent K_D • Concentration determination of active analyte in complex media (e.g. crude lysate) • Quick verification of protein expression • Contact time ~ seconds to hours so slow interactions can be measured with Static SPR • Low viscosity buffers and biomolecules 10-500 kDa are recommended² • Up to 4 independent experiments (4 channels) for P4 SPR device¹ |
| <p>Kinetic SPR</p> | <p>Purpose: Sensitive, quantitative binding characterization</p> <ul style="list-style-type: none"> • Constant flow of analyte enables ligand-analyte interaction so both strength and rate of binding are measured as Quantitative K_D • Can transfer Static SPR conditions to Kinetic SPR • Due to controlled flow, Kinetic SPR can accommodate a broader range of molecule sizes (100+ Da) • Suitable for low analyte concentrations, better real-time sensitivity |

| Data Output Differences | |
|---|--|
| <p>Static SPR</p> <ul style="list-style-type: none"> • Association without dissociation • End-point response units are analyzed • Like ELISA, Apparent K_D can be estimated at $\frac{1}{2} RU_{max}$ determined by concentration titration | <p>Kinetic SPR</p> <ul style="list-style-type: none"> • Both association and dissociation • Quantitative K_D calculated by k_{off}/k_{on} • Steeper upward curve \rightarrow faster association • Steeper downward curve \rightarrow faster dissociation |

| | Static | Kinetic |
|---|---|--|
| Use Case | Screening, end-points | Real-time binding |
| Binding parameters | Apparent K_D | k_{on} , k_{off} , K_D |
| Concentration determination | Yes | Yes |
| Analyte types & average size range | Peptides, antibodies, protein complexes (10 kDa – 500 kDa) ² | Small Molecules, peptides, antibodies, nanoparticles, vesicles (100 Da – 500+ kDa) |
| No. of experiments | Up to 4 independent experiments ¹ | 4 (Up to 2 independent experiments) |
| Amount of ligand | > 10 μ g | > 10 μ g |
| Sample injection type | Manual | Semi-automated |
| Average assay time | 5 – 15 mins | 5 – 20 mins |
| Methods used alongside SPR | Western blot, ELISA, immunoprecipitation assay, isothermal calorimetry, biolayer interferometry, fluorimetry, atomic force microscopy | |

¹Four different conditions (media, surface chemistries, ligands, etc.) in one setup using Affinité’s P4SPR device
²In some cases, small molecules and large nanoparticles can also be used - [contact us](#) for details