

## MULTICHANNEL CARTRIDGE-BASED BENCHTOP qSPR™

Leading edge bioanalytical platform for  
clinical analysis and life sciences


Affinité Instruments continues the  
portable SPR innovation with  
the world's first multichannel cartridge-  
based benchtop qSPR™

## WAIT LESS AND BIND MORE



Cartridge encases both **SPR  
sensor chip** and **microfluidic  
channels**

- ✓ Safer handling for blood, serum, urine, saliva samples
- ✓ Wifi remote access for biosafety labs
- ✓ Proprietary 4-channel design



Integrate with  
kinetics module **KNX2™**

- ✓ Real-time binding kinetics
- ✓ Insights into biological pathways and mechanisms
- ✓ Dual Channel pump